



FIGURE: PKAL-01

The patient card placement unit PKAL-01 for the exact placement of the patient card on the tray.

Product information:

- High accuracy:** The pre-printed patient cards located in the dispenser store are separated and a gripper arm removes them from the dispenser. A continuous control is carried out to ensure that only a single patient card was removed. Direct removal from a printer is also possible as an alternative.
- Positioned:** The patient card is placed on a clearly defined position on the tray.

Technical data:

Connected load:	approx. 1.3 KW
Dimensions (L x W x H):	approx. 1,300 x 750 x 1,800 mm
Required compressed air:	approx. 120 l / min.
Output:	approx. 15 cards per minute